



HIGH SECURITY BIOBANKING SOLUTION



Biobanks and biosamples storage

Biobanks collect, process, store, and distribute biological samples and their related data for a therapeutic, research or traceability purpose

Definition




OECD *OECD 2007* " A biobank is a collection of biological material and the associated data and information stored in an organized system, for a population or a large subset of a population"



Main characteristics

- A biobank can store different type of biological samples...
 - Tissues
 - Blood and derivatives components
 - DNA/RNA
 - Cell lines
 - Other biological fluids
- ... in different goals and end usages...
 - Traceability
 - Therapeutic use
 - Research on a specific pathology
 - Epidemiological research
- ...and using different storage technologies...
 - Room temperature
 - Mechanical freezer
 - Nitrogen vessels


For biosample storage, biobanks have mainly the choice between 2 types of technologies (mechanical freezer or nitrogen vessel) and 5 majors type of containers



3

Freezers mostly enable -80°C storage temperature while nitrogen enables even lower storage temperature, from $\sim -160^{\circ}\text{C}$ in vapor phase to -196°C in liquid phase

① Storage technologies




Freezers are easier to use and compatible with automated process, nitrogen offers a longer and more reliable conservation

4

Sample temperature evolution over working day

Storage technologies: Freezer

Working space temperature change over day




Front samples and back samples are not exposed to the same temperature over time

5

Autonomy of LN2 freezer after LN2 supply is stopped: 25 days

Storage technologies: LN2

















Temperature homogeneity and maintenance over time are key parameters for samples

6


Five type of consumables can be used in cryostorage , depending on the storage technology, with cryotubes, straws, high security tubes and sealed glass vials

② Storage technology relevance

Type of consumable	Freezer	Nitrogen vessel		Usage	Comments
		Vapor phase	Liquid phase		
Microplate 					Mainly short term storage, low aliquot usage flexibility, compatible with automated biobank
Cryotube 					Mass usage in freezer and LN2 vapor phase, should not be used in LN2 liquid phase
Sealed, glass vial, tube 					Glass ampule mainly used in the pharma industry, high cost of usage. High Security tube is the cost effective alternative.
Straw 					Straw are sealed and are therefore the only consumable really adapted to liquid LN2

7

Biobanking quality chain is composed of 7 major key steps. Each step is supported by SOP's. Best practices have been developed and published by international organizations (BBMRI, ISBER, OECD,...)



8

“High – quality, data – rich samples are essential for future research. But obtaining and storing these samples is not a straightforward as many researchers think”

Alison Hubel Nature / vol. 486 / 7 June 2012

Sample quality drivers

- Quality of a biological sample
 - Preservation of initial qualities, representative
 - Preservation of capacities (development,...)
 - Absence of contamination (biological, chemical,...)
 - Identification and associated data
- For future use
 - Known
 - Unknown...
- The storage technology should always be based on the future use of the biological sample: if the future use is not (fully) determined the best preservation conditions should be selected

Sample quality value

- The value of a biological sample
 - Quality of preservation (container, environment)
 - Traceability
 - Of sample origin
 - Sample process handling
 - Link with data
 - Number of aliquots
 - Minimize temperature change exposure (microplate)
 - Analytical sensitivity change over time
 - Sample exchange from lab to lab
 - Size of collection the sample belongs to
 - Mirror site storage
 - Number of research team involved (specialties & sites)
- Value can be scientific, technical, personal (patient) and/or industrial

9

“The lower storage temperature is the better”: Greater stability and long time before degradation

The 4 key quality drivers

- **Quality**
 - Only few biomolecules are preserved well at -20°C, most samples must be stored below -80°C to stop degradative enzymes and definitely at -196°C, such as live cells, at which point enzymes are thought not to be able to function at all
- **Safety**
 - Sample container must assure an **absolute security for operator** and should **avoid any cross contamination risks** between sample within the storage space
- **Traceability**
 - Essential to link sample information (Patient data, analytical data,...) to the bio specimen itself, the **sample identification is a must**
- **Efficiency**
 - To maximize sample collection size and **minimize storage cost** research should select the **right sample container**

CBS solution

- **Packaging and sample management** solutions for preservation and storage of biological samples at variable low temperatures:
 - Liquid Nitrogen phase (-196°C)
 - Liquid vapor phase nitrogen (-132°C)
 - Ultra-low temperature electrical freezer (-80°C)
- To **preserve valuable biological samples and viable cells for future use**, ensuring **maximum sample viability** by assuring the:
 - **Quality:** preserving and maintaining a sample 's initial quality
 - **Safety :** ensuring the safety of the stored product, user, other stored samples
 - **Traceability:** enhancing traceability with tamperproof linear barcoding, or 2D barcoding
 - **Efficiency:** maximizing storage space by reducing foot print



10

CBS™ High Security Straw & Co-XStraw™

- **Storage space efficiency:** up to 4-fold more aliquots vs. cryovials
- **Sample quality:**
 - Higher surface/volume ratio for better and fast freezing
 - Chemically inert and biocompatible material, no plasticizers agents including paraben, bi phenol, phthalate
 - Dnase Rnase free
 - Cryogenic grade ionomeric resin
- **Security:**
 - Full seal to eliminate risk of cross-contamination of specimens or their environment
 - Guaranteed leak-proof and shatter proof container
 - Indestructible under normal conditions of use
- **Sample identification and sample inventory management**
 - CBS™ High Security Straws are available in three different volumes: 0.3mL, 0.5mL and 2 mL (7 different colors)
 - CBS™ High Security Co-X Straw™ available on 2 different volumes: 0.3mL and 0.5mL (8 different colors)
 - Linear barcoding (EAN 128), classic straw up to 11 characters CO-X up to 18 characters
 - Goblets with 8 different colors hosting internal subdivisions
 - Canister for up to 6 goblets or racks for up to 32 goblets





11

Lineage from High Security Straw, a fully-sealed container for direct storage in liquid nitrogen in a standard dimension cryogenic vial format, with a useful volume of 1.2 mL

CBS™ High Security Tube

- **Sample quality:**
 - Takes away the barriers and safety concerns around the long term packaging of cell
 - Chemically inert and biocompatible material, no plastifying agents (paraben, bi phenol, phthalate...)
 - Dnase Rnase free
 - Cryogenic grade, ionomeric resin
- **Security:**
 - Full seal eliminates risk of cross-contamination of the specimen or its environment
 - Guaranteed leak-proof and shatter-proof container
 - Indestructible under normal conditions of use
- **Sample identification and sample inventory management**
 - Colored inserts for easy primary identification
 - 2D data matrix coded inserts
 - Compatible with cryoresistant labels
 - RFID cryoresistant chips can be placed in inserts
 - Compatible with standard racks, boxes, and canes
 - Can be used in manual or automated filling system



12



13

Filling, Sealing, Printing and Sorting final aliquot from standard sample primary tube

DIVA – CBS Automatic process

- **Walk-away design**
 - XYZ pipetting arm and consumable management
 - Automated straw distribution and conveyor
 - Straw sorting system
 - Error management and traceability
- **Autonomy**
 - 3 automated loading racks with 32 positions each
 - 8 compartment straw hopper (8 x 60 straw)
 - 2 E-Tip racks and 2 Tecan tips racks on board
 - 2 sorting drawers
- **Easy to use**
 - Touch screen user interface
 - Administrator access control
 - Production file stored on board or exported to LIMS
 - Real time production progress
- **Assured full traceability and security of produced aliquots**
 - Operator
 - Primary sample
 - Produced aliquots
- **Reduces non added value laboratory tasks**



14

Example of an upright freezer set-up with CBS racks (4 drawers) for daisy goblets, internal useful volume 680 liters: **86 016** CBS™ High Security Straw or **110 592** High Security CoX-Straw

CBS Storage solution



- CBS™ High Security Straw: **86 016 straws**
- CBS High Security CoX-Straw: **110 592 straws**



Specialty designed racks for CBS Daisy goblets


- 2 reference available
 - 022112 : 3 drawers (3 x 8); **24 goblets**
 - 021998 : 4 drawers (4 x 8); **32 goblets**
- **Straw capacity**
 - ✓ 022112:
 - CBS™ High Security Straw: **4 032 straws**
 - CBS High Security CoX-Straw: **5 184 straws**
 - ✓ 021998 :
 - CBS™ High Security Straw: **5 376 straws**
 - CBS High Security CoX-Straw: **6 912 straws**



15

Example of Liquid Nitrogen tank storage set-up with CBS canister in liquid phase, internal useful volume 1 000 liters: **211 680** CBS™ High Security Straw or **272 160** High Security CoX-Straw

CBS Storage solution



- **Daisy goblets specially designed to contain CBS High Security Straws**
 - 8 different goblet colors
 - 12 colored internal compartment, visotube
- **Straw capacity**
 - CBS™ High Security Straw: (12 x 14) **168 straws**
 - CBS High Security CoX-Straw: (12 x 18) **216 straws**
- **Canister designed to receive Daisy goblet equipped with lifter**
 - 5 different height, from 2 up to 6 goblet
 - 9 different color tab to be install on lifter

16

CBS Controlled Rate Freezer



- **Common Features**
 - Temperature range: -145°C to 45°C
 - Cooling rate: -0,1°C to -60°C per minute
 - Controller: 20 programs – 99 cooling rates – 450 segments
 - Structure 100% 304L stainless steel
- **Reliability**
 - Designed and manufactured by IMV Technologies since 1963
 - Eurotherm temperature controllers with quick release probes
- **Easy and practical use**
 - Special sample racks design (bags/tubes/glass ampules...) upon request
 - Easy to clean and disinfect
- **Traceability**
 - Automated recording of freezing curves for each cycle
 - Recordable and printable freezing process as required for GLP and SOP's
- **Security**
 - Chamber temperature alert for loading and unloading
 - Maintain -140°C at end of freezing cycle

17



18



Thank You

Michel DEPONT
Cell: +33 646 882 308
Michel.depont@cryobiosystem-imv.com

