




液態及固態乳製品檢測新科技

New Analytical Technology for Liquid and Solid Dairy Products



孫丕忠 Terence Sun  
August 26<sup>th</sup> 2020, Taoyuan, Taiwan

FOSS 2020 CONFERENCE


液態乳製品快速檢測技術

Rapid Method for The Quality of Liquid Dairy Products

FOSS 2020 CONFERENCE


FTIR 檢測技術之應用 What Can FTIR Offer Us Today?

- 自1970年代以來，中紅外線(MIR)光譜分析儀，已成為例行檢測牛奶中主要成分(脂肪、蛋白質、碳水化合物)的快速方法  
Mid-Infrared (MIR) spectroscopy has been available since the 1970's as a rapid method for routine measurements of the main constituents in milk ( fat, protein, carbohydrates )
- 傅立葉轉換式紅外光譜 ( FTIR ) 是獲取 MIR 光譜分析的最新方法  
Fourier Transform Infrared ( FTIR ) spectroscopy is the state-of-the-art method for acquiring MIR spectra analyzing



FOSS 2020 CONFERENCE


液態乳 MIR 光譜 Mid IR Spectra of Liquid Milk




- Position of the peaks → Qualitative analysis
- Intensity of the peaks → Quantitative analysis

※ Frédéric Dehareng & Clément Grelet, ICAR 2019 Conference, Technical session 7: Challenges in Creating Additional Value from Milk Analysis

FOSS 2020 CONFERENCE



※ Frédéric Dehareng & Clément Grelet, ICAR 2019 Conference, Technical session 7: Challenges in Creating Additional Value from Milk Analysis



※ Frédéric Dehareng & Clément Grelet, ICAR 2019 Conference, Technical session 7: Challenges in Creating Additional Value from Milk Analysis

ICAR 2019 會議 – FTIR 新的檢測模式  
ICAR 2019 Conference – Focus on New Tools

FOSS

Hall Panorama

<b>8:30-10:30</b>	<b>Technical Session 7 Challenges in Creating Additional Value from Milk Analysis</b> Chairpersons: Silvia Orlandini and Jere High
<b>8:30-8:50</b>	S07(T)-OP-1 Additional value of cell differentiation in the course of DHI testing Folkert Onken
<b>8:50-9:10</b>	S07(T)-OP-2 Pregnancy testing in dairy cows using a PAG test in milk samples: Different thresholds for different stages of the pregnancy Daniel M. Lefebvre
<b>9:10-9:30</b>	S07(T)-OP-3 New quality assurance challenges with recent mid-infrared models Frédéric Dehareng
<b>9:30-9:50</b>	S07(T)-OP-4 Implementation of a routine Fourier-transform infrared procedure for fatty acid analysis in milk Daniel M. Lefebvre
<b>9:50-10:10</b>	S07(T)-OP-5 Routine infrared phosphorous determination in ex-farm milk giving better insight in the phosphorous cycle on dairy farms Harrie van den Bijgaart
<b>10:10-10:30</b>	Question and Discussion

ICAR 2019 CONGRESS

MODERN UNUSUAL INNOVATIONS

ICAR 2019 會議 – FTIR 新的檢測模式  
ICAR 2019 Conference – Focus on New Tools


FOSS

<b>Electronic posters</b>	S07(T)-PP-01 „MastiMIR” - A mastitis early warning system based on MIR spectra Laura Monica Dale
	S07(T)-PP-02 Prediction of evaluated energy balance [NEL and ME] in dairy cows by milk mid-infrared [MIR] spectra Laura Monica Dale
	S07(T)-PP-03 „KetoMIR2” - Modelling of ketosis risk using vets diagnosis and MIR spectra for dairy cows in early lactation Laura Monica Dale
	S07(T)-PP-04 The use of fatty acid profiles from milk recording samples to predict body weight change of dairy cows in early lactation in commercial dairy farms Franziska Dettmann
	S07(T)-PP-05 Large scale dataset to improve and validate the prediction of lactoferrin content using milk mid-infrared spectrometry Hélène Soyeurt
	S07(T)-PP-06 A first approach to predict nitrogen efficiency of dairy cows through milk FT-MIR spectra Clément Grelet
	S07(T)-PP-07 From new milk-testing parameters to new DHI services - The view of an instrument manufacturer Daniel Schwarz

MODERN UNUSUAL INNOVATIONS

國際乳業聯盟(IDF) 公告  
Bulletin-International Dairy Federation, IDF

FOSS



IDF 447 Published 2010




IDF 504 Published 2020

SOLUTIONS FOR DAIRY ANALYSIS

FTIR 檢測模式-品質確保的新挑戰  
New Quality Assurance Challenges with FTIR Models

FOSS



同一個樣品 · 不同的光譜資料



不同的測值

不同的儀器設備

- 廠牌不同
- 同廠牌 不同機型
- 同廠牌 相同機型

同一台儀器設備

- 實驗室 溫度/濕度的影響
- 設備維修
- 設備定期保養
- 使用方式/零組件磨損

ANALYTICS BEYOND MEASURES

FTIR 檢測模式-品質確保的新挑戰  
New Quality Assurance Challenges with FTIR Models

FOSS

解決方案



光譜標準化

J. Dairy Sci. 98:2166-2169  
http://dx.doi.org/10.3168/jds.2014-8764  
© American Dairy Science Association®, 2015.

Standardization of milk mid-infrared spectra from a European dairy network

C. Grelet,<sup>1</sup> J. A. Fernández Pierna,<sup>1</sup> P. Dardenne,<sup>1</sup> V. Baeten,<sup>1</sup> and F. Dehareng<sup>2</sup>  
<sup>1</sup>Wallon Agricultural Research Center, Valoniuz of Agricultural Products Department, 24 Châteauneuf de Namur, 5030 Gembloux, Belgium

J. Dairy Sci. 100:7910-7921  
https://doi.org/10.3168/jds.2017-12720  
© American Dairy Science Association®, 2017.

Standardization of milk mid-infrared spectrometers for the transfer and use of multiple models

C. Grelet,<sup>1</sup> J. A. Fernández Pierna,<sup>1</sup> P. Dardenne,<sup>1</sup> H. Soyeur,<sup>1</sup> A. Vanlerde,<sup>1</sup> F. Colinet,<sup>1</sup> C. Bastin,<sup>1</sup> M. Gengler,<sup>1</sup> V. Baeten,<sup>1</sup> and F. Dehareng<sup>2</sup>  
<sup>1</sup>Valoniuz of Agricultural Products Department Wallon Agricultural Research Center 5030 Gembloux, Belgium  
<sup>2</sup>Agricultural Bio-Engineering and Chemistry Department, University of Liège, Gembloux Agro-Bio Tech, S135 Gembloux, Belgium  
<sup>3</sup>Wallon Breeding Association, B-5050 Ciney, Belgium




IDF 490 Published 2017

SOLUTIONS FOR DAIRY ANALYSIS



5 個生乳樣品 - Fat, BHB, Phosphorous, Nitrogen use efficiency  
6 台 MilkoScan FT+ (法國) & 7 台 MilkoScan FT 6000 (瑞士)

FOSS






MilkoScan™ FT3  
The World's First  
Intelligent Milk Analyser

### 專利光譜標準化技術-預防儀器漂移 Always Standardised-No More Instrument Drift


FOSS

- MilkoScan™ FT3 導入FOSS新取得專利之光譜自動標準化技術，此技術無需人工干預，亦不需使用FTIR標準液
- 隨著光譜自動標準化技術新專利的推出，儀器之間不再存在漂移和結果變化的情況
- MilkoScan FT3 每兩小時自動執行光譜標準化一次，消除儀器漂移，使檢測結果持續維持穩定的一致性
- 光譜自動標準化技術，不僅可以年復一年地實現檢測設備穩定的高性能，不再需要執行耗時的標準化檢查，大幅減少高成本的傳統分析校正，同時還提供了提高利潤和優化成本的機會。

現行方式光譜標準化



MilkoScan™ FT3 自動光譜標準化



ANALYTICS BEYOND MEASURES

### 自我診斷功能-確保儀器正常運作 Instrument Self-Diagnostic Ensure High Up-Time

FOSS



內建 20個 ID Chips



Warning

提供超過100個解決方案提示



ANALYTICS BEYOND MEASURES

### 檢測產品種類及項目 Full Coverage Robust Calibrations

FOSS

- Milk**  
Fat, Protein, Casein, Lactose, Low Lactose, Glucose, Galactose, Total Solids, Solids Non Fat, Density, Citric Acid, Urea, Free Fatty Acids, Titratable Acidity and Freezing Point
- Cream**  
Fat, Protein, Lactose, Total Solids and Solids Non Fat and Freezing Point
- Whey**  
Fat, Protein, Lactose, Total Solids, Solids Non Fat and Titratable Acidity

**Yoghurt & fermented application**  
Fat, Protein, Lactose, Total Solids, Solids Non Fat, Glucose, Fructose, Sucrose, Total Sugars and Lactic Acid

**Dessert & ice cream application**  
Fat, Protein, Lactose, Total Solids, Solids Non Fat, Glucose, Fructose, Sucrose, Total Sugars and Lactic Acid

**Concentrated & fortified milk**  
(Conc. Milk, Baby Milk, Infant formula, infant formula powder, Evaporated Milk, Sweetened Condensed Milk)  
Fat, Protein, Lactose, Total Solids and Solids Non Fat

**Whey concentrates & permeate**  
(WPC, WPI and concentrated permeate)  
Fat, Protein, Lactose, Total Solids, Solids Non Fat and Titratable Acidity

**Other Applications**  
Juice, Soy Sauce.....




ANALYTICS BEYOND MEASURES

- MilkoScan FT3 允許對分析樣品的粘度進行自動調整，因此可對高固體含量的WPC進行可靠的分析

### 固態及半固態乳製品快速檢測技術

Rapid Method for The Quality of Solid and Semi-solid Dairy Products

FOSS ANALYTICS BEYOND MEASURES



IDF 497 Published 2019

FOSS ANALYTICS BEYOND MEASURES

### 固態及半固態乳製品快速分析儀 FoodScan™ 2 for Solid and Semi-Solid Dairy Products Analysis

FoodScan™ 2  
A VARIANT TO SUIT YOUR NEEDS



採用 近紅外光穿透式技術 (NIR transmission) 檢測固態及半固態乳製品


例如 : butter and spreads, cheese, yoghurt, quark...等

亦可應用於肉製品檢測 (牛肉、羊肉、雞肉、魚肉...)

FOSS ANALYTICS BEYOND MEASURES

### 檢測原理 FoodScan - Operating Principle

- NIR Transmittance 400-1100 nm on a rotating sample + 1 reflectance detector, enabling color analysis
- 3 spot transmission, 3 times the volume analysed, with less rotations
- 15-60 seconds



FOSS ANALYTICS BEYOND MEASURES

1 登錄樣品資訊

2 檢測

3 得到結果

可同時測得樣品色度 (CIE L\* a\* b standard)

FOSS




REDUCE SAMPLE PREPARATION



大幅降低樣品前處理的時間



可獲得樣品更多的資訊



Accurate results every time because:

- With 3 sub-scans an even larger share of the sample is measured, which make you less dependent on the homogenisation step
- FoodScan™ 2 is a major step forward in NIR transmission technology

人工神經網路 (ANN) 技術檢量線

原物料的改變、生產製程的變化及消費者喜好的改變而產生新產品是不可避免的趨勢

以人工神經網路 (ANN) 技術建立的檢量線，具有智能性及穩固性，可確保具備處理所有問題的能力

如果所需要檢測的樣品種類或檢測項目，未涵蓋在ANN檢量線中，FOSS可提供支援協助客戶，進行客製化的服務



REDUCE SAMPLE PREPARATION

收集超過30,000筆樣品資訊建立ANN檢量線

檢量線範圍涵蓋絕大多數的固態及半固態乳製品

- CHEESE
- BUTTER
- SPREADS
- CULTURED AND FERMENTED PRODUCTS
- WHEY POWDER

優良食品安全之願景  
GoodProduct – A Strong Vision

FOSS



Tools to secure an effective and safe supply chain

REDUCE SAMPLE PREPARATION

