



# 法國IMV器材產品介紹 暨 福昌豬場精液推廣

主講人：台灣福昌集團福昌豬場 張育京



**臺灣福昌企業**  
FORMOSA FORTUNE GROUP

總公司：宜蘭市東港路22-39號

TEL: 03-9381131-5 FAX: 03-9385358 E-mail: ffgcom@ms49.hinet.net



**發昌企業有限公司**  
FA CHANG ENTERPRISE CO. LTD.

福昌豬場：宜蘭縣壯圍鄉功勞路20號

TEL: 03-9301757 FAX: 03-9304815

ISO 9001



**臺灣福昌集團**  
FORMOSA FORTUNE GROUP

## ■ 法國IMV器材暨產品介紹

## ■ BTS稀釋粉

### Short-term medium - BTS

■ BTS	QSF 1 liter (by 100)	021381
■ BTS	QSF 5 liters (by 40)	021380
■ BTS	QSF 10 liters (by 20)	021379
■ BTS	QSF 50 liters (by 4)	021378

BTS is a short-term medium formulated with high grade ingredients. Contains antibiotics.



## ■ 長效稀釋液

Highly purified BSA substitute protein complex.  
Blood protein free.

### Long-term liquid medium - SafeCell+

■ SafeCell+	100 ml	QSF 1 liter (by 10)	016181
■ SafeCell+	500 ml	QSF 5 liters (by 10)	016182
■ SafeCell+	1 000 ml	QSF 10 liters (by 10)	016183

SafeCell+ is a liquid long-term extender formulated with high grade ingredients. Contains antibiotics.



## ■ GEDIS豬冷藏精液自動分裝機

### GEDIS

NEW

#### ■ Unique all-in-one AI concept

021400

Universal, fast and convenient. Improve your staff efficiency and motivation  
No more time wasted holding the bottle or tube for minutes on end, while other sows are waiting to be bred. This innovative concept produces significant time savings during AI (up to 75%)  
The GEDIS® "auto-insemination" concept allows the technician to breed more females in a shorter amount of time while improving the quality of breeding through better heat detection and stimulation of sows.



### GEDIS® Packaging Systems

NEW

#### ■ 3 chamber machine (Serie 900) (1)

021010

Capacity: 250 / hour  
Manual placement of the Gedis® onto the machine  
Peristaltic pumps fill the Gedis®  
70 to 85 ML per dose  
Very reliable

#### ■ 2 chamber machines (Serie 400) (2)

021036

Capacity: 200 / hour



## ■ COLLECTIS 公豬精液自動採集設備

### COLLECTIS

NEW

#### Automatic boar collection system 021039

Better productivity: twice as efficient as manual collection.

Better semen quality with an insulated vagina for an optimal hygiene of collection.

Safer working conditions: no bending, no kneeling, less risks of injuries.

Easier management: lower turnover rate in your production staff with the machine.

Collectis fits all type of existing barns.



## ■ GTB1000-自動包裝機

### Automatic packaging machines

NEW

#### ■ GTB 1000: filling, sealing and labelling machine

020871

The GTB 1000 has revolutionized the perception of automatic semen packaging machines within boar studs. IMV Technologies have incorporated design features and modifications recommended by our customers.

These changes have changed the way the range of automatic GTB packaging machines looks, operates and continues to deliver high quality / high integrity semen. The GTB 1000 is compact setting the new standard in computer controlled filling systems.

The filling accuracy is significantly improved compared to all other automatic packaging machines with an improved volume tolerance of  $\pm 1g$ .

The sophisticated and sleek design is simple to use and user friendly.

Several built-in safety mechanisms will automatically stop the system in the event of a cycle failure.

The GTB 1000 will fill / seal and label in excess of 1,000 bags per hour. One person is capable of not only running the machine but also of performing additional tasks while the GTB are being packaged.

Runs with 6-bar compressed air - delivered without compressor





# 恒溫水浴槽

## Water-bath

■ 30 l water bath with stand for baby flasks and tubes (4)

007227

■ Pyrex baby flasks with caps (by 24)

020177



## 精子濃度測定儀

### Accuread: precision counts

The only portable solution for semen analysis

Accurate

Species specific adaptation

Easy to use

Lightweight unit


Low maintenance

Easy to clean

Battery/100 to 240 V



## 自動分析精子型態儀



## 顯微鏡


### Microscope binocular objectives 10x - 20x - 40x

- Binocular microscope YS100  
with heated stage and histology [1] 016745
- Binocular microscope YS100  
with heated stage without histology 016744



### Microscope trinocular objectives 10x - 20x - 40x

- Trinocular microscope E50i  
with heated stage and histology [2] 019172
- Objective 100X achromat plan for microscope E50i 019171
- 36 cm monitor with color camera 016513
- Slides 76x26 mm (by 10x50) 005795
- Cover slides 18x18 mm (by 1,000) [3] 006557
- Mini pipette for semen sampling (by 10x100) 007283



## 冷凍精液保存相關設備



- Mini-Digitcool for 24 racks [7] 007263

- Freezing rack for 58 straws 0.5 ml [8] 007118


- Distribution block for 58 straws 0.5 ml 007092

- 60 l tank with flexible + casters + gauge  
for Mini- & Micro-Digitcool 020086

- ← ■ 120 l tank with flexible + casters + gauge [3] 007223

- ← ■ 200 l tank with flexible + casters + gauge [4] 007222





# 謝謝指導

台灣福昌集團

總公司:宜蘭市東港路22-39號

電話:03-9381131

精液訂購:03-9301538

郵箱: [ffgcom@ms49.hinet.net](mailto:ffgcom@ms49.hinet.net)