



Features

iSperm is designed for vet, simplifying the process of artificial insemination and production control. It provides a user-friendly app to accomplish sperm analyzing precisely and efficiently.

iSperm Animal Sperm Analyzer



Function

Motility: Total Motility

Analysis Time: 10 sec video analysis

Progressive motility: LIN/STR/VSL/VAP/VCL

Concentration: Total Concentration, Dilution Doses

Sample: Fresh semen, Diluted semen, Frozen semen

Supported Species: Swine / Poultry / Canine / Caprine / Bovine / Equine



Web Site QR Code

Contact

Taiwan:
P: +886-2-23685358
Email: service@aidmics.com

China:
P: +86-756-8933120/+86-1581260869
Email: chiang@aidmics.com

Global:
P: +886-2-23685358
Email: service@aidmics.com